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New Human Derived Product Advanced for Cell Culture and Stem Cell therapy

AventaCell BioMedical Co., Ltd. ("AventaCell") is among the world leaders devoted to developing novel human-derived products for use in cell culture and tissue regeneration. Helios Bioscience is AventaCell's product family brand for their products to be used in cell culture and tissue regeneration. AventaCell technologies used in the Helios Bioscience line offer new human-derived solutions for use in translational research of cell and tissue-based therapies to meet the need for animal serum-free cell expansion and production. The demand for safe, efficient and cost-effective cell expansion and production is rapidly increasing with the growth in cell therapy and regenerative medicine research and clinical development. Helios Bioscience products are designed to support expansion and production of a broad range of cells including mesenchymal stem cells and multiple immunocell lines. AventaCell is committed to providing animal serum-free products to accelerate the research, development and commercialization of safe and efficacious cell and tissue-based therapeutics.

Non-Xenogenic Supplement Animal Serum Free

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Virus Category	RNA Enveloped	RNA Non-Enveloped	DNA Enveloped	DNA Non-Enveloped
Model for Virus	HCV, HIV	HAV	CMV, EBV, HBV	B19
Family	BVDV	Reo3	HSV1	MMV
Genome	Flavi	Reo	Herpes	Parvo
Size (nm)	ssRNA	dsRNA	dsDNA	ssDNA
Resistance	40 - 60	60 - 80	120 - 200	18 - 24
UltraGRO-PURE GI (Log10 reduction)	Low	Med-High	Medium	Very high
	<b>&gt; 5.42</b>	<b>&gt; 4.40</b>	<b>&gt; 4.51</b>	<b>4.55</b>

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

Fig1:

Gamma irradiation has low impact on MSCs expansion potency of UltraGRO™-PURE (UG-P). Media with 5% (v/v) gamma irradiated UltraGRO™ supplements retained comparable growth kinetics and yields.

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Table2: Immuno-phenotypical characterization of human MSCs. Human MSCs derived from adipose tissue (AD), umbilical cord matrix (UC), bone marrow (BM) cultured in UltraGRO™-PURE GI for 5 passages display characteristic expression of MSC surface markers.

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
Fig.2: Human adipose tissue derived MSCs retain tri-lineage differentiation capability after cultured in UltraGRO™-PURE GI supplemented medium for three passages.

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

Fig.3: MSCs retain immunomodulation potency

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